



	biotema tıbbi cihazlar ltd. şti.		
	Dumlupınar mah. Doruk Sk. ISTANBUL 216 ikinci etap C BLOK Kat 9 No 70 34720 KADIKÖY		
	biotemaVATS Instruments		
Product Name	Photo	Model	Description
Atrial Needle - Holding Forceps (Head with one - side open)		RXQCZ004	shaft diameter 8mm, curved head, total length 340mm
		RXQCZ005	shaft diameter 8mm, straight head, total length 340mm
		RXQCZ006	shaft diameter 6mm, curved head, total length 340mm
		RXQCZ007	shaft diameter 6mm, straight head, total length 340mm
Thoracic Tissue Forceps (Head with one-side open)		RXQZZ004	shaft diameter 8mm, head width 10mm,lymph forceps, total length 280mm
		RXQZZ011	shaft diameter 6mm, head width 10mm, lymph forceps, total length 280mm
		RXQZZ017	shaft diameter 6mm, head width 6mm, lymph forceps, total length 280mm
Thoracic Tissue Forceps (Head with one-side open)		RXQZZ009	shaft diameter 8mm, head width 10mm, oval forceps(with teeth), total length 280mm
Thoracic Tissue Forceps (Head with one-side open)			shaft diameter 6mm, head width 10mm, oval forceps(with teeth), total length 280mm
		RXQZZ010-430	head with one - side open,shaft diameter 6mm, head width 10mm, oval forceps(with teeth), total length 430mm
		RXQZZ012	shaft diameter 6mm, head width 8mm,oval forceps(with teeth), total length 280mm
Thoracic Tissue Forceps (Head with one-side open)		RXQZZ013	shaft diameter 6mm, head width 6mm, oval forceps(with teeth), total length 280mm
		RXQZZ014	shaft diameter 6mm, head width 10mm, oval forceps(without teeth), total length 280mm
Thoracic Tissue Forceps (Head with one-side open)		RXQZZ015	shaft diameter 6mm, head width 10mm, C forceps(without teeth), total length 280mm
Thoracic Tissue Forceps (Head with one-side open)		RXQZZ016	shaft diameter 8mm, head width 12mm, oval forceps(with teeth), total length 280mm
		RXQZZ018	shaft diameter 6mm, head width 8mm, oval forceps(without teeth), total length 280mm
Ligature forceps(Alis forceps) (Head with one-side open)		RXQJZ005	shaft diameter 8mm, head width 5mm, total length 340mm
Thoracic Hemostatic Forceps (Head with one-side open)		RXQZX026	shaft diameter 6mm, curved 15° , total length 340mm
Thoracic Hemostatic Forceps (Head with one-side open)		RXQZX024	shaft diameter 5.5mm, curved 30° , total length 330mm
		RXQZX019	shaft diameter 6mm, curved 30° , total length 340mm
		RXQZX005	shaft diameter 8mm, curved 30° , total length 340mm
Thoracic Hemostatic Forceps (Head with one-side open)		RXQZX020	shaft diameter 6mm, curved 60° , total length 340mm, right hand
		RXQZX020-340	shaft diameter 6mm, curved 60° , total length 340mm, left hand
		RXQZX020-430	DeBakey head width 5mm, shaft diameter 7mm, total length 430mm
		RXQZX012	shaft diameter 8mm, curved 60° , total length 340mm
Thoracic Hemostatic Forceps (Head with one-side open)		RXQZX021	shaft diameter 6mm, curved 90° ,right angle, total length 340mm
		RXQZX016	shaft diameter 8mm, curved 90° ,right angle, total length 340mm
Thoracic Hemostatic Forceps (Head with one-side open)		RXQZX022	shaft diameter 6mm, curved 90° ,small right angle, total length 340mm
		RXQZX018	shaft diameter 8mm, curved 90° ,small right angle, total length 340mm
Thoracic Hemostatic Forceps (Head with one-side open)		RXQZX023	shaft diameter 6mm, curved 90° ,big right angle, total length 340mm
		RXQZX014	shaft diameter 8mm, curved 90° ,big right angle, total length 340mm
Thoracic Hemostatic Forceps (Head with one-side open)		RXQZX027	shaft diameter 6mm, curved 70° ,slant right angle, total length 340mm
Thoracic Forceps (Head with one-side open)		RXQZX025	shaft diameter 6mm,straight head, total length 340 mm
Blocking Forceps		RXQZX028	shaft diameter 8mm, curved head, double row non-traumatic teeth,total length 340 mm
Cardiac Surgery Forceps (Head with one-side open)		RXQJD003	shaft diameter 6mm, curved head, length 320mm
		RXQFY007	shaft diameter 7mm, head width 8mm, total length 340mm
		RXQFY012	shaft diameter 5.5mm, head width 5mm, total length 340mm

Triangle Lung Grasping Forceps (Head with one-side open)		RXQFY008	shaft diameter 7mm, head width 5mm, total length 330mm, right hand
		RXQFY008-330	shaft diameter 7mm, head width 5mm, total length 330mm, left hand
		RXQFY008-430	Snake head width 5mm ,shaft diameter 7mm ,total length 430mm
		RXQFY009	shaft diameter 7mm, head width 10mm, total length 310mm
Auricle clamp (Atrial appendage forceps)		RXQZX030	shaft diameter 6mm, head length 6mm, tip of the head 3X3, bending length 105°, height 44mm, non-traumatic teeth, total length 330mm
Tissue Clamps(Thoracic Tissue Forceps)		RXQZZ016-430	Oval forceps head width 12mm, shaft diameter 8mm, total length 430mm
Thoracic Tissue Forceps (Esophageal Anastomosis Clamp)		RXQZZ019	shaft diameter 10mm, head length 50mm, tip of the head 5X10, height 28mm, non-traumatic teeth, total length 340mm
Purse Stitching Clamp		RXQZZ020	double joint, 5-9 teeth, total length 360mm
Knot Pusher		HJH01095	shaft diameter 5mm, Y- type shallow, total length 330mm
		HJH01095V	shaft diameter 5mm, Y- type deep, total length 330mm
Knot Pusher		RXWDJ001	shaft diameter 5mm, C type, total length 360mm
		RXWDJ002	shaft diameter 5mm, V type, total length 360mm
Electric Hook		FQD901A	straight, L-shaped head, shaft diameter 5mm, working length 330mm
		FQD901B	straight, L-shaped head, one curved, shaft diameter 5mm, working length 330mm
		FQD901C	straight, L-shaped head, one curved, shaft diameter 5mm, working length 330mm
		FQD901E	straight, 7-shaped head, shaft diameter 5mm, working length 330mm
Suction with Electric Coagulation		FQC721C	shaft diameter 5mm ,curved hollow coagulation, total length 390mm
Flushing Suction Tube (Two Ways Shank)		RXWXY001	(1) + (6) + (3), (1) + (7) + (3), (1) + (8) + (3), 5mm shaft diameter, total length 360mm, metal head, complete suction, you can choose the items from the former three combination
		RXWXY003	(1) + (6) + (4), (1) + (7) + (4), (1) + (8) + (4), 5mm shaft diameter, total length 360mm, metal head, complete suction, you can choose the items from the former three combination
		RXWXY003-450	One curved head(1 big hole, around 4 small holes), metal handle, shaft diameter 5mm, blunt tip, should be smooth, without knurling, no frosted, total length 450mm
		RXWXY004	(1) + (6) + (5), (1) + (7) + (5), (1) + (8) + (5), 5mm shaft diameter, total length 360mm, metal head, complete suction, you can choose the items from the former three combination
Dennis Suction Instrument		RXWXY003-455	5 mm shaft, curved blunt tip, working length: 12-1/3" (32.0 cm), total length: 17-3/4" (455 mm)
Flushing Suction Tube (Two Ways Shank)		RXWXY001-01	(2) + (6) + (3), (2) + (7) + (3), (2) + (8) + (3), 5mm shaft diameter, total length 360mm, plastic head, complete suction, you can choose the items from the former three combination
		RXWXY003-01	(2) + (6) + (4), (2) + (7) + (4), (2) + (8) + (4), 5mm shaft diameter, total length 360mm, plastic head, complete suction, you can choose the items from the former three combination
		RXWXY004-01	(2) + (6) + (5), (2) + (7) + (5), (2) + (8) + (5), 5mm shaft diameter, total length 360mm, plastic head, complete suction, you can choose the items from the former three combination
Flushing Suction Tube		RXWCX001(HJH01092B)	shaft diameter 5mm ,straight head, total length 360mm

(Three Ways Shank)		RXWCX002	shaft diameter 10mm ,straight head, total length 360mm
Flushing Suction Tube (Four Ways Shank)		RXWCX003	shaft diameter 5mm ,straight head, total length 360mm
Trocars		RXWCC001	diameter 10mm, straight, with thread
		RXWCC002	diameter 10mm, straight, without thread
		RXWCC003	diameter 5mm, straight, with thread
		RXWCC004	diameter 5mm, straight, without thread
Foerster Lung Grasping Clamp		RXQFY009-424	7 mm shaft, ring handle, curved right (left-hand style), 5mm oval ring jaws, working length: 7-1/2" (19.0 cm), total length: 11-1/2" (424 mm)
Diego Gonzalez-Rivas		RXQFY009-419	7 mm tapered shaft, ring handle, mixer-style 1×2 DeBakey jaws, with ratchet, working length: 12-1/2" (32.0 cm), total length: 16-1/2" (419 mm)
Node Grasper		RXQZZ016-330	5 mm shaft, ring handle, curved left, 11 mm oblong ring jaws, with ratchet, total length 330mm
Node Grasper		RXQZZ016-424	5 mm tapered to 7 mm shaft, ring handle, curved left, 11 mm oblong ring jaws, with ratchet, total length 424mm
Dennis Rib Cutter		RXQZZ016-10-16	10 mm blade height, 16 mm cutting width, straight shaft, working length 290mm, total length 340mm
Dennis Rib Cutter		RXQZZ016-15-25	15 mm blade height, 25 mm cutting width, straight shaft, working length 290mm, total length 340mm
Dennis Rib Cutter		RXQZZ016-20-25	20 mm blade height, 25 mm cutting width, straight shaft, working length 290mm, total length 340mm